



For a better environment

Prescription for total Bio - Protection
& Sterile Environment



Professional We can

BADGER-D

The disinfection power of colloidal silver & Hydrogen Peroxide

Eco Friendly for Hospital &
Pharma Application

Aerial Fumigant, Surface and Water Disinfectant Solution

Composition:

Colloidal Silver IP 0.01% w/v

Hydrogen peroxide IP 10% w/v

Activity:

Broad spectrum: Bactericidal, Fungicidal, Virucidal, Sporicidal and Tuberculocidal

Features & Benefits:

Broad - spectrum - kills all pathogens including spores and resistant strains like MRSA and VRE Just takes 1 hour to sterilize OT. Can be used to mop OT surfaces. Excellent surface and material compatibility.

Made from Entirely Food safe ingredients and Biodegradable. Most effective against Bio film.

Applications:

MEDICAL:

In hospitals and laboratories for aerial fumigation in critical areas such as OT, maternity units.

FOR AERIAL FUMIGATION USE G SAKTHI FOGGER FOR BETTER RESULTS

Use 1:20 BADGER i.e., 200ml BADGER + 800ml DM water

INDUSTRY:

In pharmaceutical, food and dairy industry as a surface disinfectant for aerial fumigation in manufacturing, R&D and QA/QC areas.

Usage Directions :

AERIAL FUMIGATION

Use 1:20 BADGER i.e., 200ml BADGER + 800ml DM water

SURFACE DISINFECTION:

Use 1:5 BADGER i.e., 50ml BADGER + 950ml DM water

WATER SANITATION:

Use 100 ppm i.e., 100ml of BADGER in for 1000 Litre of water water with contact time of 8 hours

Use 500 ppm i.e., 500ml of BADGER in for 1000 Litre of water water with contact time of 6 hours

Use 1000 ppm i.e., 1 Ltr of BADGER in for 1000 Litre of water water with contact time of 2 hours



Other Technical information, available on request :

Basil Hygiecare

95, Balan Nagar, Peelamedu, Coimbatore - 641 004

Phone : +91 422 2593600 | Mobile : +91 9655771984

E-Mail : info@basilhygiecare.com | www.basilhygiecare.com



For a better environment